


www.elva-1.com sales@elva-1.com

USER OPERATION AND MAINTENANCE MANUAL

D-band one channel Interferometer Part No. MMI-140R2



2013

Table of Contents.

- Introduction.
 1.1 General Description.
- Specifications.
 - 2.1 Electrical.
 - 2.2 Mechanical.
 - 2.3 Block-diagram of interferometer
 - 2.4 Principle of operation
- Installation and Functional Tests.
 3.1 Assembly procedure.
 3.2 Test and tuning procedure.
- 4. Measurements.
 - 4.1 Put into operation
 - 4.2 Measurements.
 - 4.3 View and Save data.
 - 4.4 Software tools.
- 5. Test report.
 - 5.1 Output Power
 - 5.2 Operation Frequency.
 - 5.3 AGC signal.
 - 5.4 Accuracy of measurement.
 - 5.5 Test measurement of real time signal

1. INTRODUCTION.

This instruction manual contains information on installation and operation of the D-band Interferometer.

1.1 General Description.

D-band interferometer is intended for measuring the line-averaged density of a plasma along the path through which the mm-wave beam is passed, through phase shifts in the propagated beam.

Base principle of operation is an effect of change of phase speed of electromagnetic waves of a millimeter range in plasma depending on density.

2.1 Electrical Specifications.


 RF Frequency LO Frequency Frequency Stability Time Phase Analyze Time collection Phase error measurement Data rate of real time density Real time density interface Real time signal rate Output RF Power Max attenuation Tx-Rx IQ Detector Frequency RF to IF Gain Conversion Loss of Mixer Noise Figure Gain of Horn Waveguide Flange Data interface 	140.0 GHz; 138.0 GHz; 10 ppm; 5 μs; 1000 ms; 0.07 radian; 10kHz voltage, +0.1 +5.0V 15.4 Rad/V 100 mW; > 90dB 100 kHz; 70 dB with AGC system; 9 dB; 12 dB; 30 dB; WR-06 UG-387/U-M TCP/IP Ethernet; TTL;
18. Flange 19. Data interface	UG-387/U-M TCP/IP Ethernet;
20. External trigger21. Input impedance of trigger input22. Waveguide transition23. AC Power24. Operating temperature range:	50 Ohm; WR-06 to WR-42; 220 VAC (external modules) +10°C+ 50°C;

2.2 Mechanical Specifications.

- 1. Transmitter
- 2. Receiver

300x205x70 mm; 250x235x70 mm;





2. 4 Principle of operation.

The interferometer operates at fixed frequency. Due to all built-in oscillators are phaselocked by one reference Cristal oscillator 100 MHz, there is a possibility to control phase shift of mm-wave signal is passed through plasma.

Mm-wave oscillators (transmitter and LO for receiver) are built on one principle. The transmitter consists of ultra low phase noise DRO 7.0 GHz, direction coupler 1:10, doubler and multiplier x10 with integrated high power amplifier at 14 GHz. The direction coupler provides 7.0 GHz signal for receiver. The transmitter generates output frequency at 140.00 GHz, 7.0 GHz x20.

The mm-wave receiver consists of VCO 6.9 GHz, direction coupler 1:10, doubler, multiplier x10 with integrated high power amplifier at 14 GHz and balance mixer. For stabilization IF signal and phase-locking LO frequency, the signal 7.0 GHz from transmitter is used. After down conversion the signal from transmitter we have IF frequency 2.0 GHz. Then this signal is amplified by LNA and comes to power amplifier with built-in AGC system. AGC system allows to keep stable amplitude of IF signal in 30-70 dB attenuation range in mm-wave channel. Response time of AGC system is 2.5 microsec. Then stabilized IF signal is applied to IQ mixer, where as LO oscillator synthesizer 2.00010 GHz is used. The synthesizer 2.000010 GHz is also phase-locked by reference Cristal oscillator 100 MHz. IQ mixer provides two IF signals 100 kHz (sin and cos), which are digitized by built-in 12 bits ADCs. The controller converts digital data in Ethernet pockets and transfers them via TC/IP interface to PC station.

PC station collects data from the interferometer and calculates absolute phase changing during discharge of plasma. The software operates the following way: as 100 kHz signals and trigger of ADC converters are locked by one reference, we can control only pass zero points of IF signals without any additional reference signal, i.e. relation between period of IF signal and quantity of digitized points inside of this period is constant in anytime (without plasma). Such way we get absolute rate Pi / (digitized points) φ o and can control phase shift of signal during discharge of plasma. As we have two 100 kHz signals from IQ detector (sin and cos) total pass zero points is four times more and finally time phase analyze is 2.5 microsec. The software calculates time different between close pass zero points and compare it with absolute rate. Then sum different between actual and absolute rate:

$$\sum_{N}^{0} \left(\varphi(n) - \varphi(n-1) \right)$$

Where,


n - number of current pass zero point from the beginning of discharge. N - total quantity pass zero points during of discharge

After collection data and calculation total drift of phase the software recalculate phase in absolute changing of density of plasma with rate $0.4143*10^{(18)}$ / m³ (operation frequency 140GHz, diameter of plasma 0.4m)

3. INSTALLATION AND FUNCTIONAL TEST.

3. 1. Assembly procedure.

The interferometer can be installed directly near windows of TAKAMAK or removed from it by extension WR-42 waveguides, as it presented on picture below.



Wave guide layout sceme for Interferometer


Waveguide's set of the interferometer consists of the following components:

- Waveguide transition WR-06/ UG-387/U-M to WR-42 / UG-595/U





- Horns WR-06/ UG-387/U-M





After setting waveguides the following connection according block-diagram of the interferometer should be done:





1. 7 GHz reference signal from the transmitter (1) to receiver by coaxial cable with SMA connectors. The system is completed 2 cables, one is spare.



2. Power supplies via applied cables with DB-9 connectors. Power supplies are equal and can be used any for transmitter and receiver.



3. Receiver to PC station by TP5 cable with RJ-45 connector.



4. Special cable for control IF signals 100 kHz and AGC signal.




- 5. BNC input for external trigger TTL pulse
- 6. SMA output 'Real time density', 0 ...+5V

Connect all BNC connectors any oscilloscope for control 100 kHz IF signals and AGC.

When all connections are done, complete PC station with standard monitor, keyboard and mouse and switch ON all power supplies and PC.

3. 2. Test and tuning procedure.


Before installation the system on TOKAMAK it can be tasted on table. Assemble the system as disrobed above and place Tx/Rx modules as presented on picture below:



Then run program 'Interferometer' (placed on desktop), the main menu of program appears:

Interferometer	×
Measure	
View	
Test	
Quit	

For test IF signals 100kHz and AGC signal test mode should be run, press button 'Test' and real time oscilloscope window with 100 kHz IQ signals appears



For test trigger operation Measure mode should be run, press button 'Measure' and new sub-menu is appearing:


Connection	
IP address 192.168.10.50	UDP port 63
Output Discharge no. 129 Folder C:\Users\800378\Deskt	op\data\Browse
	Measure Clos

then press button Measure in sub menu for waiting a trigger pulse. Enter discharge No and select folder for saving data.


Apply trigger pulse (or make short cut) to BNC connector (5). After collection of data come back to main menu (press Close) and then select View mode.

Папка:	🌗 data		•	← 🗈 💣 💷 ▼	
C:	Имя	*		Дата изменения	Тип
Недавние места Рабочий стол	00129_201	30912_173640.ifr		12.09.2013 17:36	Файл "IFF
Библиотеки					
(Щ) Компьютер					
	•	III			۲
Сеть	Имя файла:	00129_20130912_173640		•	Открыть
	Тип файлов:	Raw data (*.fr)		-	Отмена

Select file with set before discharge No and open it, after that measured result is appeared. Picture below demonstrates correct operation of the system.



After passing 'table test' the system can be installed on TOKAMAK. Assemble and adjust complete system and before starting to work the user has to test and tuning horns (if it is needed). For that the user should switch on system in Test mode and check form of IF signals and level of AGC system. To control sinusoidal wave form of IQ signals and AGC signal tune horns such way to get minimum available level of AGC signal. Below dependence of AGC signal vs attenuation in waveguide channel is presented. For reference, total attenuation in assembling system with installed horns in open area should be about 60 dB, that corresponds to about +2.95V output of AGC signal.




After tuning procedure the user should test trigger operation of system as it was described above.

4. MEASUREMENTS.

4. 1. Put into operation.

1. Connect cable with trigger signal from TOKAMAK control system to opto-isolated BNC connector (5) on the receiver. The trigger pulse should be TTL standard and has duration no less then 1 ms. The interferometer starts to collect data with rise front of the trigger pulse.

2. Run software and put into menu 'Measure'



3. Enter discharge serial No (1) and select folder (2) for saving data and then press Measure (3) for regime 'waiting a trigger'

IF Measurer	X J
Connection IP address 192.168.10.50	UDP port 63
Output Discharge ng. 129 Folder C:\Users\800378\Desktop\data\	Browse 3 Measure Close

4. The system is ready for measurement.

4.2. Measurements.

After start up procedure the system is waiting for trigger. When trigger pulse is appeared software automatically save raw data into selected folder, add '+1' to discharge serial No and then comes back into regime 'waiting for trigger'. To stop circle of collection data, press button Close.

UDP port	63	
8		Browse
	1	
	ke.	Measure

5. Software generates file name with real time indication:

00129_20130912_173640.ifr

where,

discharge No_year_moth_day_hour_min_sec

The file structure is:

•••••


Where, N – quantity of measurements No of measurement _ ADC 16 bits code I channel _ ADC 16 bits code Q channel


4. 3. View and Save data.

To view saved data come to main menu or run one more Interferometer software and then press button View. Selection file window is appeared.

ipen file		and the second		X
Папка:	퉬 data	-	🗢 🗈 💣 📰 🕈	
Ca	Имя	*	Дата изменения	Тип
Недавние места	00129_201	30912_173640.ifr	12.09.2013 17:36	Файл "IFF
Рабочий стол				
Г Библиотеки				
(Щ) Компьютер				
6	٠	III		•
Сеть	Имя файла:	00129_20130912_173640	•	Открыть
	Тип файлов:	Raw data (*.ifr)	-	Отмена


Select file and press open, when software recalculates raw data and draw result.






User can switch On/Off presented data to use flags into window No1. Below is example with measured data only.


The result of measurement can be presented as a phase shift in radians (Y-axis is graduated in radians) or electron density (Y-axis is graduated in $*10^{(18)} / m^{3}$). To see result in radians in window No2 '1.000' should be entered (picture above) and '0.4143', as presented below






There is possibility to install markers at measured dada. For that press marker list, activate flag 'add markers' (window No3) and then select desired point by using PC mouse tool.


To save data as image in bmp format, press 'Export image', enter file mane and select folder for saving.



To save final data in txt format press 'Export data', enter file mane and select folder for saving. Software stores data in file 'xxxx.ifd'




4. 4. Software tools.



5. FUNCTIONAL TESTS.

5. 1. Output power.

Test setup in presented below.





Measured data: Pout = 104.4 mW

5.2. Operation frequency.


Test setup in presented below.



Measured data: Fout = 140,004271 GHz


5.3. AGC signal.

Test setup in presented below.






Test data is presented on plat below.



5.4. Accuracy measurement.


Test setup in presented below





Measured data: Accuracy measurement is < 0.07 radians

5.5 Test Measurement of real time signal

Test setup in presented below









Principle of operation – apply voltage to electromagnetic solenoid and move reflector. Coil spring returns reflector in initial position. Micrometer check actual shift of reflector. Trigger oscillator provides pulse for apply voltage to solenoid and run the system.


Measured data

Micrometer: 8.8mm, that corresponds to 8.8 / 2.14 (wave length for 140GHz) *6.28 =25.8 rad Real time signal on the oscilloscope:





The system measured 25.97 rad



D-band Interferometer